




CS ELITE

WALL MOUNTED CHEST STAND SYSTEM

The XR Fusion Medical CS Elite system is designed for facilities that specialize in Pulmonology and Occupational Medicine. Advanced features like ergonomically correct release handle, electric locks and remote collimator light make positioning fast and efficient.

- 30-50kW generators
- Extensive vertical travel – 16”- 69”
- Accommodates a full range of upright exams
- Pair with an optional mobile table to expand imaging capabilities
- DR upgradeable - ask about our digital upgrade solutions



www.brownsmedicalimaging.com

Brown's Medical Imaging
14315 C Circle Omaha, NE 68144
1 (800) 701-XRAY (9729)

Prestige Medical Imaging
289 N Plank Rd, Suite 4
Newburgh, NY 12550
(845) 561-6947



XR-LITR045-0523